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
EN

Before using, please read this manual completely and save it for future use. On our website you will find the user manual for the Streamcheck App.

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1.1. Intended Purpose and intended Users

Streamcheck is a medical self-testing device for males aged 18 and older, enabling urinary health assessments at home. It supports

long-term monitoring and evaluation of urinary function through repeated measurements. The app generates structured PDF reports that can be shared with healthcare professionals. The device itself does not diagnose or interpret data but provides reliable measurements to support clinical decision-making. To ensure proper functioning of the device, it must be used only under the following environmental conditions: 0 °C – 30 °C; 10 % - 90 % RH; 700 hPa – 1060 hPa.

Core Functions

- Urine Flow Monitoring: Records flow patterns over time with a maximum deviation of $\pm 10\%$ to assess bladder function.
- Irregularity Detection & Guidance: The app automatically analyzes results. If abnormalities are found (e.g., low Qmax or pathological biomarkers – color-coded green, orange, red, gray), it provides clear instructions, such as repeating the test or consulting a doctor.
- Data Output: After each test, the app generates a structured PDF report containing the flow curve, biomarker results, timestamps, and traceability data. This report is for clinical interpretation only and contains no diagnostic statements or automated evaluations.

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EN 1.2. Contraindications

The Streamcheck device must not be used by users:

- with severe cognitive impairments that prevent understanding and following the instructions for use, unless used under supervision.
- with pronounced motor impairments (e.g. paresis, severe tremor) that make safe and effective handling of the device impossible.
- with known allergies to materials used in the device housing that come into direct skin contact (e.g. ABS or MABS plastics).

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2. Safety instructions and safety related warnings

All information in this part is provided to instruct you to use the appliance in a correct and safe way so that you can avoid injuries to yourself as well as damage to property. If any questions arise please contact the manufacturer.

- Always read the instructions for use and the user manual in full before use.
- Adherence to the IFU is mandatory to ensure proper and effective operation.
- If you have any health concerns regarding the measurement results of the Streamcheck, consult a urologist. This device is used to record medical data. No diagnoses are made.
- The use of the device does not replace medical consultation and treatment. Therefore, always consult your doctor first if you experience any type of pain or illness.
- Observe the general hygiene measures when using Streamcheck.
- Keep children away from the packaging material (risk of suffocation).
- Do not charge it unattended. Remove the charger when the Streamcheck device is fully charged to maximise battery life and minimise the risk of malfunction.
- Do not immerse the Streamcheck device in water. Protect the device from strong impact and vibration.
- Ensure that the device is operated at a sufficient distance from sources of electromagnetic radiation, such as WLAN routers, mobile

phones, microwave ovens and other radio devices. A minimum distance of 1.5 metres should be maintained to ensure the proper functioning of the Streamcheck device

- If the Streamcheck device has been dropped or is visibly damaged, it shall no longer be used. If in doubt, contact customer service or your distributor.
- Any serious incident that has occurred in relation to the device should be reported to Streamcheck GmbH and the competent authority of the EU Member State in which the device was acquired.
- The Streamcheck may only be operated with the corresponding cups. The use of third-party cups will result in incorrect measurements and may damage the device.
- Only use the cup with urine test strips once. The urine test strips only work the first time they are used and cannot provide any results thereafter.
- Only use the device: on people, for the purpose for which it was developed and in the manner specified in these instructions for use.
- Under no circumstances may you open or repair the device, as otherwise its proper functioning can no longer be guaranteed. Failure to do so will invalidate the warranty.
- This device must not be modified without the manufacturer's permission. If the device is modified, the manufacturer can no longer ensure its safe use.
- The appliance is maintenance-free.
- There is a risk of measurement distortion due to incorrect positioning (e.g. penis resting on the cup ring).
- There is a risk of malfunctioning due to external influences or environmental conditions such as magnetic fields, external electrical and electromagnetic effects, electrostatic discharge, radiation associated with diagnostic or therapeutic procedures, pressure,

humidity, or temperature. Always store your device in the provided case.

- Measurements below a total urine volume of 150 ml may not provide valid uroflowmetry results; the app will inform the user accordingly.
- In case of symptoms suggesting an active urinary tract infection (e.g. fever, abdominal pain), users are advised to consult a physician before using the device.

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







EN 3. Start-up and operation

1. Download the Streamcheck app (Apple-Store / Google Play).
2. Insert the cup into the Streamcheck device so that the mark in the cup is aligned with the mark in the Streamcheck ring. The QR code will then be on the right side.
3. Now press the cup under the holding clamps. To do so, apply slight counter-pressure to the ring holding the cup.
4. Start the Streamcheck by touching the Streamcheck logo on the device for about 2 seconds. Open the 'Streamcheck' app on your smartphone.
5. Click "Add measurement" in the app and the app will start pairing the device. The process is displayed in the app. The completed pairing is indicated by a solid blue LED light (see 4. LED displays, item 10).

6. The QR code on the cup is automatically read in (detected by the camera) and stored in the app for single use. The LED light will display a white running light as confirmation (see 4. LED displays, item 5).
7. Urinate into the cup. Make sure that nothing is placed on the cup and that the cup itself is not resting on anything that could affect the measurement.
8. Empty the cup by tipping it out and place it on its side on a firm surface. When placing the cup, avoid direct sunlight.
9. Once the measurement is complete, the results can be viewed in the smartphone app.
10. To remove the cup from the Streamcheck device, squeeze the cup in your hand until the circumference is smaller than the ring and slowly slide it out of the device. Finally, dispose of the used cup as advised in chapter "05. Cleaning and Disposal", and clean the cup ring with the provided wipes.
11. After a successful analysis of the measurement and data transfer, the device will be automatically switched off.

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EN 5. Charging

- Use the charger provided to charge the Streamcheck. Connect it to a CE-certified power adapter with a USB-C port (5 V | 1000 mA).
- The charger can remain permanently inside the storage case. This way, the Streamcheck charges while in the case, without the need

to remove the charger. At the same time, you ensure that the device is always correctly positioned.

- The charging process starts when the LED display lights up orange. When fully charged, the LED display will switch off.
- The device's battery level is also displayed in the app.
- The full charging process takes about 1.5 hours.
- The charging status can also be seen from the LED display (see 04. LED-Display).

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EN 6. Cleaning and Disposal

Cleaning

The outer casing of the Streamcheck and the glass protection of the camera must be carefully cleaned and disinfected after each use. Use only the cleaning wipes provided with the case. Do not immerse the device in liquid.

Disinfection using UV radiation is not permitted.

Waste disposal

The product must be disposed of in accordance with the waste management regulations. Electronic components must be disposed of in accordance with the regulations for waste electrical and electronic equipment.

Package disposal

The plastic packaging around the cups should be disposed of in the yellow recycling bin.

Cup disposal

The used cups must be disposed of with the residual waste.

Disposal of used appliances

In the interest of environmental protection, the device must not be disposed of with household waste at the end of its lifespan. Disposal can be done through appropriate collection centres in your country.

Dispose of the device in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, contact the municipal authority responsible for waste disposal.

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EN 7. Troubleshooting

Hereby are shown some potential issues and how to solve them.

Issue	Solution
Device returns to pairing mode after use. The device has restarted.	Wait approx. 2 minutes, the device will automatically switch to standby.
The device does not switch on when the Streamcheck logo is touched.	Charge the device.
The device does not pair with the smartphone.	Close the application, re-activate Bluetooth and restart the application.
The cup cannot be recognised.	Check that the cup is correctly aligned in the Streamcheck holder.

Issue	Solution
The application cannot connect to the device.	Switch Bluetooth off and on again on your mobile phone. Make sure that the Streamcheck device is charged. See 4. LED Display point 2

If any issues persist or you have any other inquiries, please contact the customer support.

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EN 8. Declaration of Conformity

This product meets the requirements of the EU Medical Device Regulation MDR 2017/745. The undersigned, acting on behalf of the manufacturer, hereby declares that the product described above meets the relevant requirements of the EU Medical Device Regulation (MDR) 2017/745.

This product has been designed, manufactured and tested in accordance with the following EU directives and standards (EU Medical Device Regulation MDR 2017/745).

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Number: 0494

Conformity assessment procedure: CB-Zertifikat EMV

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All measurements in the App are stated in SI units. Volume is shown in ml (with a tolerance of +/- 5%) and Uroflow is shown in ml/s (with a tolerance of +/- 10%)

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3x: Protection against solid objects over 2.5mm in diameter; x2: Protection against vertically dripping water when tilting the device to 15 degrees.

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